

# CURIE

## Versatility that Delivers



### Enlightening your creations

CURIE is equipped with a sophisticated UV-light source management system that enables users to do more in less time and space. Ideal for dental professionals who take bio-compatibility seriously and care about patient safety.



#### Efficient Curing Space

Perfect curing chamber space for Dental Clinics and Labs with rigorous workflows.



#### Future-proof

Equipped with an SD card reader for updating material settings and post-curing times.



#### Accuracy and consistency

Preloaded optimal curing times database available at the turn of a dial and touch of a button.



#### Premium Performance

Smart UV-light source management system for maximized energy efficiency and the ideal light strength to prevent any deformation.

### TECHNICAL SPECIFICATIONS

<b>Dimensions</b>	20 (W) x 22 (D) x 16 (H) cm 7.9 (W) x 8.7 (D) x 6.3 (H) in
<b>Chamber Dimensions</b>	14 (W) x 14 (D) x 8 (H) cm 5.5 (W) x 5.5 (D) x 3.1 (H) in
<b>Weight</b>	Approx. 3 kg Approx. 6.6 lbs
<b>Power Consumption</b>	75 W
<b>UV Wavelength</b>	365, 385, 405 nm
<b>Power Input</b>	100-240VAC, 50-60Hz
<b>Curing Time</b>	1 sec - 20 min

### PACKAGE INCLUDES

- 1 x CURIE UV Oven
- 1 x Glass Plate (Inside the machine)
- 1 x Power Cable
- 1 x SD Card
- 1 x Card Reader



Scan the QR Code to  
Schedule a Virtual Demo Today!



**ORDER NOW!**