

# SO | Multi | Units | A-dec® Orthodontic



# Set your practice apart with A-dec.

If you have a vision to set your orthodontic practice apart, precision matters. Aesthetics matter. Craftsmanship matters. Not just in the treatment you provide, but the equipment you use.

A-dec orthodontic solutions are thoughtfully designed, down to the smallest detail. Efficient functionality and proven performance come together with thoughtful enhancements and refined touches. All in a complete, customizable solution that enhances the treatment experience for you, and your patients.

This is orthodontic equipment made like no other. Because it's made only by A-dec.



# DESIGNED FOR THE WAY YOU WORK

A-dec orthodontic solutions work together seamlessly for a smooth and easy workflow—day after day, patient after patient. Because when you move efficiently, so does your practice.

- Patient Chairs
- Specialty Workstation & Delivery
- Dental Stools
- LED Lights
- Sterilization Center
- Mechanical Room Equipment

And when it comes to your well-being, you can feel the difference with A-dec. Our proven, ergonomic designs help you do more than just work efficiently—they enhance your comfort, support good form, and help you maintain a healthy posture so you can work longer and more comfortably.







# MAKE THE MOST OF YOUR SPACE

For many orthodontic practices, an open bay floor plan is the optimal way to keep your team moving efficiently from patient to patient throughout the day. Placing rows of chairs and workstations side by side uses less square footage than a traditional layout—and supports a smooth flow of patients in and out of treatment.

A-dec orthodontic solutions feature small footprints and compact designs that help you optimize an open bay setup within your available square footage. Choosing A-dec helps you make the most of your valuable floor space—without having to compromise on equipment quality, comfort or performance.





# ENGINEERED FOR COMFORT

A comfortable patient is a happy patient, so A-dec designs every dental chair with that in mind.

Sleek, armless designs make it easy to get in and out of the chair. Pressure-mapped, contoured surfaces—and a virtual pivot that moves with the patient's body—create a relaxed experience. A dual-articulating, gliding headrest and ample shoulder support keep patients resting comfortably.

And with a thin, flexible backrest that lets you roll in close and maintain a healthy posture, A-dec dental chairs are just as comfortable for you as they are for your patients. Choose from sewn or formed upholstery, and customize the color to complement your practice design.



## A-dec 317 Specialty Chair

- Chair recline switch with automatic back up
- Fixable height setting from 19.5" to 16.5"(495 mm to 419 mm)
- Optional mounting for A-dec LED lights
- Supports up to 350 lb (159 kg)\*
- 5 year warranty



## A-dec 311 Dental Chair

- Mechanical lift for smooth, flexible positioning
- Base height range from 31.5" to 13.75" (800 mm to 349 mm)
- Swivels 30° each side of center
- Optional mounting for A-dec LED lights
- Supports up to 400 lbs (181 kg)\*\*
- 5 year warranty



\*Load rating based on UL requirement of 4X (1400 lb/636 kg) overload static testing.

\*\*Load rating based on UL requirement of 4X (1600 lb/725 kg) overload static testing.

# ORGANIZATION COMES NATURALLY

When there's a place for everything, it's easy to get into an uninterrupted treatment flow.

The A-dec Specialty Workstation accommodates everything you need—equipment, technology, and materials—all in one place. The versatile design offers ample drawer and countertop space, including room for essential integrated equipment. The articulating delivery arm has a 10” swing-out range to keep handpieces within easy reach, with no twisting or straining.

Choose from two cabinet sizes, customizable with the full palette of A-dec Inspire® solid surfaces, quartz, and laminates. Then add an LED light, monitor mount, ICV capability, or other cabinet options to suit your practice needs.

## A-dec 397 Specialty Workstation

- Solid surface countertop
- Five full-extension, soft-close drawers
- Computer storage with fan
- Subbase with easy access to utilities

## A-dec 341sp Specialty Delivery

- 5-position delivery configuration (includes saliva ejector and high volume ejector)
- Optional touchpad
- Wet/dry foot control
- Quick-disconnect air and water

## 28" Specialty Workstation

shown with added options:  
quartz countertop, 2 glove dispensers, trash drop, light/monitor mount brackets (light and monitor sold separately), ICV cutout (ICV sold separately), standard touchpad, power supply connector kit.

## 24" Specialty Workstation



# COMPLETE SOLUTIONS FOR YOUR PRACTICE

A-dec offers reliable, creative solutions for orthodontic practices of every size.

## LED Lights

Touchless, flexible, and crystal clear, our award-winning A-dec LED dental lights are brighter than halogen lights, but use 80-90% less power. All with a life expectancy of 20 years. Choose from two models, both with the flexibility to mount to the patient chair or the Specialty Workstation.

## Sterilization Centers

Everything about A-dec sterilization centers is designed to intuitively and visually follow established CDC guidelines for instrument processing protocols, while streamlining the flow of materials. And with the ability to configure the furniture, you can design the workflow to fit your exact space and practice needs.


## Dental Stools

A-dec's dynamic seating system helps your spine remain in a healthy S-shape—the natural curve of your body—to help reduce disc pressure and muscle strain. Strategically engineered to conform and move with you, our easily adjustable stools support both comfort and improved performance.

## Mechanical Room Equipment

Designed for the dental environment, A-dec QuietCore+™ oil-free compressors and dry vacuums deliver powerful performance with remarkably quiet operation. On-demand technology channels air and suction exactly where needed, while the compact, stackable designs are ideal for smaller mechanical room spaces.





View in 3D: Rotate, zoom and pan to see all angles of your customized operator.

[a-dec.com/dream](http://a-dec.com/dream)



# THE DIFFERENCE IS IN THE DETAILS

With A-dec, quality is built in from the start. All product design, engineering, manufacturing, and assembly is done in-house—we even produce the majority of our own parts from raw materials. That attention to detail is just one reason why A-dec equipment performs so reliably and lasts so long.

Simple, intuitive designs lower the risk of failure and reduce the need for maintenance. Controlling the quality of each component down to the last bolt ensures that A-dec equipment will always meet our exacting standards. And yours.

And because A-dec equipment is tested for 20 years of life, and backed with a five-year warranty, you can choose it with confidence.




# ALIGNED FOR YOUR SUCCESS

Your A-dec ownership experience starts with installation and training—and that's just the beginning. Certified A-dec product experts are available 12 hours a day, five days a week to answer your questions and help you make the most of your equipment. Together with your A-dec Territory Manager, they provide all the support you'll need. Not just for the life of our products, but the life of your practice.

Want to learn more about A-dec orthodontic solutions?  
Get connected with your local A-dec Territory Manager:








# Reliable Creative Solutions

FOR DENTAL PROFESSIONALS WORLDWIDE

 @A-decDental

 @A-decDental

 @A-dec Inc.

 @AdecDental

 @A-dec Inc.




## More than a warranty. It's a commitment.

At A-dec, our focus is on making dentistry better—which includes superior customer service along with superior dental products. Your A-dec equipment is backed by a five-year warranty, and a lifetime of customer support. We'll be here today, tomorrow, and well into the future. That's our commitment to you.

Any colors and products shown may vary and are subject to change without notice. Contact your authorized A-dec dealer for the most current product information. All trademarks in this document are subject to the trademark rights of A-dec, its affiliates, and their collective suppliers in the United States and foreign countries.

©2024 A-dec Inc. All rights reserved. 85.1075.00/E-only/10-24/Rev B2



A-dec Headquarters  
2601 Crestview Drive  
Newberg, Oregon 97132 USA  
P: 1.800.547.1883 within USA/CAN  
P: +1.503.538.7478 outside USA/CAN

**a-dec.com**

Other Locations  
A-dec Nashville Experience Center  
A-dec Australia  
A-dec United Kingdom

