Scaling

Air Polishing

Cavitron® Tabletop System Comparison





Cavitron' Plus Cavitron' 300 Series





Cavitron[®] Prophy Jet[®]



Cavitron[®] Jet Plus

	Portable & Compact	Improved Comfort & Efficiency	The Uultimate Scaling Experience	Air Polishing	Combo Scaling & Air Polishing
Scaling and Air Polishing Capabilities					✓
Prophy Mode Auto Cycles				✓	✓
Wireless Tap-On Technology Foot Pedal		✓	✓	✓	✓
Rinse Function		✓	✓	✓	✓
Purge Mode		✓	✓	✓	✓
Pivoting Handpiece Cable	✓	✓	✓	✓	✓
Water Control on Handpiece	✓	✓	✓	✓	✓
Boost Mode	✓	✓	✓		✓
Intelligent Feedback Control (SPS Technology)	✓	✓	✓		✓
Steri-Mate/Jet-Mate Autoclavable Detachable Handpiece	✓	✓		✓	✓
Steri-Mate 360 Autoclavable Detachable Swivel Handpiece			✓		
Touch Screen Interface			✓		
Memory Presets			✓		
Featherweight Detachable Handpiece Cord			✓		
Water Filter Indication Light			✓		
Optional Reservoir Pump/Bottle System (NOTE: SKU# 90146 is sold separately and is only compatible with Cavitron Select SPS)	SKU# 90146				
Warranty	2 Years	2 Years	3 Years	2 Years	2 Years
Dentsply Sirona Part Number	81302	8184001	8270310	8195001	8187501



Package Options

Cavitron®

Tabletop Package Comparison Select™ SPS™

	Select™ SPS™ Package	Plus Package	300 Series Package	Prophy Jet Package	Jet Plus Package
Warranty	3 Years	3 Years	4 Years	3 Years	3 Years
Cavitron Care Maintenance Plan	Not Available	2 Years	3 Years	2 Years	3 Years
Total Handpieces	1 Steri-Mate Handpiece	1 Steri-Mate Handpiece	3 Steri-Mate 360 Handpieces	1 Jet-Mate Handpiece	3 Jet-Mate Handpieces
Cavitron Inserts (NOTE: (1) SlimLine 10S FitGrip included in package, Redeem 3 of your choice unpon registration. Not applicable to the Prophy Jet Package)	4 Cavitron 30K Ultrasonic Inserts	4 Cavitron 30K Ultrasonic Inserts	4 Cavitron 30K Ultrasonic Inserts	2 Cavitron Air Polishing Inserts	4 Cavitron 30K Ultrasonic Inserts + 2 Cavitron Air Polishing Inserts
Dentsply Sirona Part Number	8161424	8161425	8161429	8161427	8161428



Cavitron Care Maintenance Plan

Available exclusively with Cavitron Packages!

Cavitron Care will send you a maintenance reminder annually to schedule service on your Cavitron for the following Items:

General Operation and Function:

- Evaluation and verification of all clinical functions
- Evaluation and repair of malfunctions
- Installation of recommended factory upgrades

	~ .
Electrical	System:

- Calibration of ultrasonic circuits
- Test, verify and repair all electrical functions as needed
 Inspection and repair of
- electrical cable connections

Unit Model
Cavitron 300 Series (G310)

Cavitron Jet Plus (G137) Cavitron Plus (G136)

Cavitron Prophy Jet (G138)

Cavitron Select SPS (G124)

• Inspection, repair/replacement

water flow with replacement

• Evaluation and/or replacement

components found to be leaking

of water control valves

Testing and calibration of

of clogged or damaged

of water supply line filter

· Repair or replacement of any

components

Water Systems:

Scheduled Service Included

3 years
3 years
2 years
2 years
Not Available

Air System: (Jet Plus & Prophy Jet Systems)

- Evaluation and/or replacement of control valves
- Replacement of internal tubing and filters
- Testing and calibration of air flow rate with replacement of needed components
- Evaluation and/or replacement of main air supply line filter
- Repair or replacement of any components found to be leaking air

Cavitron[®] Ultrasonic Inserts

Each Cavitron Package includes (1) Slimline 10S FitGrip insert Choose (3) additional inserts from the options below!

Slimline[®]



10 - Left Curve Low to Medium Power 30K 80396 30K Fitgrip^{*} 82006



10 - Right Curve Low to Medium Power 30K 80397 30K Fitgrip* 82007



10 - Straight Low to Medium Power 30K 80395 30K Fitgrip" 82005



1000 Low to High Power 30K 81570 30K Fitgrip* 82008



Powerline™

Low to High Power 30K 80797 30K Fitgrip' 82001



10 Low to High Power 30K 80294 30K Fitgrip^{*} 82002



100 Low to High Power 30K 80798 30K Fitgrip' 82003



1000 Low to High Power 30K 80799 30K Fitgrip* 82004

Simplified registration and redemption at: cavitronsystems.com

When selecting inserts, consider the following:

Inserts are uniquely made to handle different levels of deposit. Consider the force needed to remove the deposits.



Thicker inserts with larger tips are designed for heavier levels of calculus.

Thinner inserts with slim, adaptable tips are designed for the removal of light to moderate deposits or biofilm.

Treatment site

Select inserts that are specially designed for the type of tooth anatomy you will be working on.



Straight inserts are designed for flat tooth anatomy.

Curved inserts are designed to give you the best results in furcations and concavities.

Thinsert*



Thinsert Low to High Power 30K 81551 30K Fitgrip' 82009



© 2020 Dentsply Sirona Cavitron° is a registered trademark of DENTSPLY International and/or its subsidiaries. 90696-1 Rev. US CCO 0319

Type of deposit Inserts are unique