Planmeca PlanMill®60 S

For advanced CAD/CAM solutions.

The most sophisticated chairside milling unit on the market.

The versatile Planmeca PlanMill 60 S unit is a dry or wet mill capable of handling discs and blocks. Equipped with a high-speed spindle and an automatic bur changer for 10 tools, the powerful milling unit has been designed specifically for high-quality restorations and prosthetics.

PRODUCT FEATURES

- 5-axis chairside unit with lab capabilities for high production
- · Choice of wet or dry milling
- · Mill discs or blocks
- Broad range of material options
- · C-clamp adapter for more detailed milling
- Automatic bur changer for 10 tools
- Open architecture for easy integration

The Planmeca PlanMill 60 S produces dentures, full arch restorations, night guards, surgical guides, and up to 25 crowns in a single disc.



Wet milling



Dry milling



Automatic tool changer for up to 10 burs



Planmeca PlanMill® 60 S













Mill blocks or discs with a broad range of materials.