




---

# Primescan<sup>®</sup> 2

Powered by DS Core<sup>®</sup>

Enter a new era of digital patient care





# Primescan<sup>®</sup> 2

## The first cloud-native intraoral scanning solution

Our direct-to-cloud, wireless intraoral scanner that stands for simplicity and versatility, helps you grow your practice, and deliver excellent patient care.



### Direct to cloud scanning

Primescan 2 is beyond wireless, it scans in the cloud and is not attached to a specific computer. DS Core<sup>1</sup>, Dentsply Sirona's cloud platform, enables you to scan with any device with an internet connection, like your operator PC, dental chair computer or tablet, allowing you to initiate a scan any time and anywhere.



### Increase scanner availability

As soon as the scan is complete, the scanner can be moved to the next patient, while consultation, treatment planning or even chairside design can be done on any computer, independently from the scanner device. This allows for a fast hand-over of the scanner within the practice and increases the availability of your device.



### One scanner for numerous treatments

Primescan 2 delivers outstanding accuracy and reliable scan quality no matter the indication or material. Thanks to its broad use, it can open a world of possibilities for your practice offering – such as restorative, implantology, orthodontic and preventive treatments.



### Proven technology<sup>3</sup> with improved hardware design

Speed, accuracy and ease of use, brought to you in an ergonomic shape: The new design includes a sleeve with a lean profile for even more patient comfort, while the scanner hardware is easy to handle for you and your team.



### Caries detection

Early diagnosis is essential for preserving oral health. Using fluorescence and near-infrared technologies, Primescan 2 can assist with caries identification<sup>2</sup>, while visualization with DS Core helps your patients to understand their treatment need.



Are you ready to enter a new era of digital patient care?

Learn more about Primescan 2 and schedule a demonstration

<sup>1</sup> Primescan 2 requires an active, paid subscription to DS Core

<sup>2</sup> Please refer to your Dentsply Sirona representative to confirm availability in your region

<sup>3</sup> Primescan's precision and accuracy have been proven in several studies. For more detail, please refer to the Dentsply Sirona Clinical Study Overview

